

One-Piece Cuff Sizes and Ranges

		Range (cm)	
Size		Minimum	Maximum
Non-FlexiPort™	Neo 1 (Size 1)	3.3	5.6
	Neo 2 (Size 2)	4.2	7.1
	Neo 3 (Size 3)	5.4	9.1
	Neo 4 (Size 4)	6.9	11.7
	Neo 5 (Size 5)	8.9	15

FlexiPort™	Small Infant (Size 6)	7	10
	Infant (Size 7)	9	13
	Small Child (Size 8)	12	16
	Child (Size 9)	15	21
	Small Adult (Size 10)	20	26
	Adult (Size 11)	25	34
	Adult Long (Size 11L)	25	34
	Large Adult (Size 12)	32	43
	Large Adult Long (Size 12L)	32	43
Thigh (Size 13)	40	55	

