


Welch Allyn Cuff Part Number Logic

Example



Part Number Breakdown

STYLE		DASH	SIZE		DASH	NUMBER OF TUBES		No Space	CONNECTOR/FITTINGS	
Description	1ST SECTION		Description	2ND SECTION		Description	3RD SECTION		Description	4TH SECTION
Soft Disposable	SOFT	-	Sm. Infant*	06	-	1 Tube	1		Screw	SC
Vinyl Disposable****	VINYL	-	Infant	07	-	2 Tube	2		Bayonet*	HP
Durable	REUSE	-	Small Child	08	-	Blank Cuff			Tri-Purpose**	TP
			Child	09	-				Locking	MQ
			Small Adult	10	-				Locking (Mated)***	MF
			Adult	11	-				Inflation System	BV
			Adult Long	11L	-					
			Large Adult	12	-					
			Large Adult Long	12L	-					
			Thigh	13	-					

*Only Available in One-Tube Configurations

**Not Available in Vinyl Disposable Cuffs. Also, add 'E' for longer tube version (ex. TPE)

***Male and Female Connectors (Two-Tube Only)

****Not available outside USA

