

TRIMLINE

Connector and Adapter Guide

Connectors

REF 1360
Package of 10 threaded plastic female screw connectors for 5/32" tubing (GE/Critikon - Dinamap monitors).
WA Equivalent 5082-164



REF 1360N
Package of 10 threaded plastic female screw connectors for 1/8" tubing (GE/Critikon Neonatal Monitors).



REF 1376
Package of 10 threaded male swivel collar connectors for 5/32" tubing, plastic.
WA Equivalent 5082-174



REF 66018
Package of 10 threaded male swivel collar connectors for 1/8" ID tubing, plastic.
WA Equivalent 5082-175



REF 66024
Package of 10 threaded male swivel collar connectors for 3/16" ID tubing, plastic.



REF 1355
Package of 1 each, threaded male cable connector for 5/32" tubing, metal (plated brass).
WA Equivalent 5082-161



REF 1361
Package of 10 quick connect male bayonet Rectus for 5/32" tubing, plastic
(Philips/HP, Spacelabs, Datascope, Siemens)
WA Equivalent 5082-184



Connectors

REF 1319

Package of 10 plastic bayonet type female quick connectors for 5/32" I.D. tubing, Hewlett Packard style.



REF 1365

Package of 1 each, metal bayonet type female quick connector for 1/8" I.D. tubing, Hewlett Packard style.

WA Equivalent 5082-177



REF 1362


Package of 10 female locking luers, lug style slip for 5/32" tubing, plastic (Datascope, Colin, Spacelabs).



REF 1367

Package of 10 plastic male locking luers, thread style with barbed end for 5/32" tubing.

WA Equivalent 5082-165



REF 1369

Package of 10 plastic male subminiature swivel locking (quick coupling) connectors with barbed end for 5/32" I.D. tubing.

WA Equivalent 5082-178



REF 1370

Package of 10 plastic female subminiatures (quick coupling) (no swivel) for 5/32" I.D. tubing (GE/Marquette, Protocol).

WA Equivalent 5082-182



REF 1371

Package of 10 slip type male luers, plastic, for 5/32" I.D. tubing.

WA Equivalent 5082-166



Connectors

REF 1372

Package of 10 male to male connectors, plastic, for 5/32" tubing.

WA Equivalent 5082-170



REF 1373

Package of 10 neonatal lug style female slip luers, plastic, for 3/32" I.D. tubing.



REF 1377

Package of 10 plastic 5/32" to 1/8" reduction fitting connectors.

WA Equivalent 5082-171



REF 1378

Package of 10 standard male slip luers with barbed end (metal) for 5/32" I.D. tubing.

WA Equivalent 5082-167



REF 1391

Package of 10 standard female slip luers with barbed end (plastic) for 5/32" tubing.



REF 1379

Package of 10 plastic ANSI/HIMA slip type male connectors, for use with 3/32" I.D. tubing, neonatal.

WA Equivalent 5082-176



REF 1380

Package of 10 plastic ANSI/HIMA slip type male connectors, for use with 1/8" I.D. tubing.




REF 2863


Package of 10 plastic male micro luers with barb, for use with 3/32" I.D. tubing, neonatal.




Connectors are combined with Tubing to make Adapters



'Y' Tube: in 1 Tube to 2 Tube Orientation



'Y' Tube: in 2 Tube to 1 Tube Orientation



'H' Tube: 2 Tube to 2 Tube




'I' Tube: 1 Tube to 1 Tube




Coiled Extension Tubing

Adapters are used to convert a Device to accept a Cuff it otherwise would not accept.



MONITOR – 2 Tube Mated
Subminiature Connectors
(Datex/Ohmeda)



ADAPTER – Double Mated
Subminiature to Single
Female Bayonet (Rectus)



CUFF – 1 Tube Bayonet
(Rectus) Male

IMPORTANT NOTES:

- **SAP and Price List Descriptions describe what fittings are on the adapter, NOT what fittings plug into the adapter.**
- **The first fitting(s) in the description describes the “MONITOR” side.**
- **The second fitting(s) in the description describes the “CUFF” side.**
- **Monitor end of adapters have a tether attached used to secure adapter to monitor tubing.**

EXAMPLE FOR ADAPTER ABOVE:

Catalog Number: 9333
Price List Description: Adapter, Double Mated Subminiature to Single Female Bayonet (Rectus), 1/EA
SAP Description: ADAPTER,DBL MF to SGL FHP

Understanding SAP Abbreviations

Catalog Number: 9333
 Price List Description: Adapter, Double Mated Subminiature to Single Female Bayonet (Rectus), 1/EA
 SAP Description: ADAPTER,DBL MF to SGL FHP
40 Character limit forces use of abbreviations in SAP descriptions

SAP Abbreviation	Definition	SAP Abbreviation	Definition
BARB	Hose Barb (Hose Barbs are always Male)	GGE	Bainbridge Pocket Aneroid Gauge
BT	Bare Tube (Bare Tubes are always Female)	HP	Male Bayonet (Rectus) (Quick Coupler)
BV	Bulb & Valve	MF	Mated Subminiature
BVBTC	Bulb & Valve, Bare Tube and Clamp	ML	Male Luer Slip
BVFL	Bulb & Valve and Female Luer Slip	MLC	Male Luer Slip and Clamp
BVHPC	Bulb & Valve, Male Bayonet (Rectus) Connector and Clamp	MML	Male Micro Luer
BVSCC	Bulb & Valve, Female Screw Connector and Clamp	MQ	Female Subminiature
DBL	Double	MSC	Male Screw (Threaded) Connector
FHP	Female Bayonet (Rectus) (Quick Coupler)	NK	Male Subminiature (Nihon Kohden)
FL	Female Luer Slip	SC	Female Screw (Threaded) Connector
FLBT	Female Luer Slip & Bare Tube	SCC	Female Screw (Threaded) Connector and Clamp
GAUGE	Bainbridge Pocket Aneroid Gauge	SGL	Single

Adapters for connecting cuffs with: 1 Bare Tube (1/8" ID)

Catalog Number	PRICE LIST DESCRIPTION
9001	Adapter, Double Female Threaded Connector to Single 1/8" Hose Barb , 1/EA
9011	Adapter, Double Mated Subminiature to Single 1/8" Hose Barb , 1/EA
9021	Adapter, Double Female Subminiature to Single 1/8" Hose Barb , 1/EA
9041	Adapter, Single Female Subminiature to Single 1/8" Hose Barb , 1/EA
9051	Adapter, Single Bayonet (Rectus) to Single 1/8" Hose Barb , 1/EA
TL9211	Adapter, Single Bare Tube to Single 1/8" Hose Barb , 1/EA




Adapters for connecting cuffs with: 1 Bare Tube (5/32" ID)



Catalog Number	PRICE LIST DESCRIPTION
9282	Adapter, Double Female Threaded Connector to Single 5/32 Hose Barb , 1/EA

Adapters for connecting cuffs with: 1 Tube with 5/32" Hose Barb




Catalog Number

PRICE LIST DESCRIPTION

9324

Adapter, Double Male Luer Slip to Single Bare Tube, 1/EA

Adapters for connecting cuffs with: 1 Tube Female Bayonet (Rectus)



Catalog Number


PRICE LIST DESCRIPTION

66054

Adapter, Double Female Luer Lock to Single Male Bayonet (Rectus), 10/PK

Adapters for connecting cuffs with: 1 Tube Male Bayonet (Rectus)

Catalog Number	PRICE LIST DESCRIPTION
9314	Adapter, Double Female Threaded Connector to Single Female Bayonet (Rectus), 1/EA
9321	Adapter, Double Female Subminiature to Single Female Bayonet (Rectus), 1/EA
9327	Adapter, Double 5/32" Hose Barb to Single Female Bayonet (Rectus), 1/EA
9328	Adapter, Single Male Subminiature to Single Female Bayonet (Rectus), 1/EA
9329	Adapter, Single Female Threaded Connector to Single Female Bayonet (Rectus), 1/EA
9333	Adapter, Double Mated Subminiature to Single Female Bayonet (Rectus), 1/EA
9335	Adapter, Double Male Luer Slip to Single Female Bayonet (Rectus), 1/EA
9338	Adapter, Double Male Subminiature to Single Female Bayonet (Rectus), 1/EA
9343	Adapter, Double Bayonet (Rectus) to Single Female Bayonet (Rectus), 1/EA
9345	Adapter, Single Female Subminiature to Single Female Bayonet (Rectus), 1/EA
9346	Adapter, Double 1/8" Hose Barb to Single Female Bayonet (Rectus), 1/EA
9222	Adapter, Single Female Subminiature to Single Female Bayonet (Rectus), Metal, 1/EA
9223	Adapter, Double Female Subminiature to Single Female Bayonet (Rectus), Metal, 1/EA
9243	Adapter, Single Male Luer Lock to Single Female Bayonet (Rectus), Metal, 1/EA
9286	Adapter, Single Female Threaded Connector to Single Female Bayonet (Rectus), Metal, 1/EA
9287	Adapter, Double Male Luer Slip to Single Female Bayonet (Rectus), Metal, 1/EA
9288	Adapter, Double Female Threaded Connector to Single Female Bayonet (Rectus), Metal, 1/EA
9297	Adapter, Double Mated Subminiature to Single Female Bayonet (Rectus), Metal, 1/EA
9344	Adapter, Single 5/32" Hose Barb to Single Female Bayonet (Rectus), Metal, 1/EA
9358	Adapter, Single Bare Tube to Single Female Bayonet (Rectus), Metal, 1/EA
TL9225	Adapter, Single Female Luer Lock to Single Female Bayonet (Rectus), Metal, 1/EA



Adapters for connecting cuffs with: 1 Tube Male Slip Luer



Catalog Number	PRICE LIST DESCRIPTION
9263	Adapter, Single Bayonet (Rectus) to Single Female Luer Slip , 1/EA
9266	Adapter, Single Female Luer Lock to Single Female Luer Slip , 1/EA
9337	Adapter, Single Bare Tube to Single Female Luer Slip, 1/EA
9341	Adapter, Double Male Luer Slip to Single Female Luer Slip, 1/EA
9352	Adapter, Double Female Threaded Connector to Single Female Luer Slip, 1/EA

Adapters for connecting cuffs with: 1 Tube Male Subminiature



Catalog Number	PRICE LIST DESCRIPTION
9295	Adapter, Single Bayonet (Rectus) to Single Female Subminiature , 1/EA
9319	Adapter, Double Female Threaded Connector to Single Female Subminiature, 1/EA

Adapters for connecting cuffs with: 1 Tube Female Luer Slip



Catalog Number	PRICE LIST DESCRIPTION
9003	Adapter, Double Female Threaded Connector to Single Male Luer Slip, 1/EA
9013	Adapter, Double Mated Subminiature to Single Male Luer Slip, 1/EA
9023	Adapter, Double Female Subminiature to Single Male Luer Slip, 1/EA
9043	Adapter, Single Female Subminiature to Single Male Slip Luer, 1/EA
9053	Adapter, Single Bayonet (Rectus) to Single Male Slip Luer, 1/EA
9073	Adapter, Single Female Threaded Connector to Single Male Slip Luer, 1/EA
9143	Adapter, Single Male Subminiature to Single Male Slip Luer, 1/EA

Adapters for connecting cuffs with: 1 Tube Female Screw (Threaded) Connector



Catalog Number

PRICE LIST DESCRIPTION

66027	Adapter, Single Female Luer Lock to Single Male Threaded Connector, 10/PK
66033	Adapter, Single Male Luer Lock to Single Male Threaded Connector, 10/PK
66038	Adapter, Single Bare Tube to Single Male Threaded Connector, 10/PK
66041	Adapter, Single Bayonet (Rectus) to Single Male Threaded Connector, 10/PK
66158	Adapter, Single Female Subminiature to Single Male Threaded Connector, 10/PK
9290	Adapter, Single Bayonet (Rectus) to Single Male Threaded Connector, 1/EA
9300	Adapter, Single Male Subminiature to Single Male Threaded Connector, 1/EA
9307	Adapter, Double Female Threaded Connector to Single Male Threaded Connector, 1/EA

Adapters for connecting cuffs with: 1 Tube Male Neonatal Luer Slip



Catalog Number

9323

PRICE LIST DESCRIPTION

Adapter, Single Male Micro Luer to Single Female Neonatal Luer Slip, 1/EA

Adapters for connecting cuffs with: 2 Bare Tubes (1/8" ID)




WA PART- NUMBER	PRICE LIST DESCRIPTION
9002	Adapter, Double Female Threaded Connector to Double 1/8" Hose Barb , 1/EA
9022	Adapter, Double Female Subminiature to Double 1/8" Hose Barb , 1/EA
9042	Adapter, Single Female Subminiature to Double 1/8" Hose Barb , 1/EA
9052	Adapter, Single Bayonet (Rectus) to Double 1/8" Hose Barb , 1/EA
9072	Adapter, Single Female Threaded Connector to Double 1/8" Hose Barb , 1/EA
9162	Adapter, Double 1/8" Hose Barb to Double 1/8" Hose Barb , 1/EA
9322	Adapter, Single Male Subminiature to Double 1/8" Hose Barb, 1/EA
9012	Adapter, Double Mated Subminiatures ro Double 1/8" Hose Barb

Adapters for connecting cuffs with: 2 Bare Tubes (5/32" ID)



WA PART- NUMBER	PRICE LIST DESCRIPTION
9283	Adapter, Double Male Luer Slip to Double 5/32" Hose Barb, 1/EA
9348	Adapter, Single 1/8" Hose Barb to Double 5/32" Hose Barb, 1/EA
9356	Adapter, Double Mated Subminiature to Double 5/16" Hose Barb, 1/EA

Adapters for connecting cuffs with: 2 Tube, 5/32" Hose Barbs




**WA PART-
NUMBER**

PRICE LIST DESCRIPTION

9234

Adapter, Double Female Subminiature to Double Bare Tube, 1/EA

Adapters for connecting cuffs with: 2 Tube, Male Bayonet (Rectus)



WA PART- NUMBER	PRICE LIST DESCRIPTION
9349	Adapter, Single Female Subminiature to Double Female Bayonet (Rectus), 1/EA
9350	Adapter, Single Bayonet (Rectus) to Double Female Bayonet (Rectus), 1/EA
9351	Adapter, Single Female Threaded Connector to Double Female Bayonet (Rectus), 1/EA
9353	Adapter, Single Female Luer Lock to Double Female Bayonet (Rectus), 1/EA
9355	Adapter, Single 5/32" Hose Barb to Double Female Bayonet (Rectus), 1/EA
9306	Adapter, Single Bayonet (Rectus) to Double Female Bayonet (Rectus), Metal, 1/EA


Adapters for connecting cuffs with: 2 Tube, Male Luer Lock

**WA PART-
NUMBER**


66042

PRICE LIST DESCRIPTION

Adapter, Single Male Luer Slip to Double Female Luer Lock, 10/PK



Adapters for connecting cuffs with: 2 Tube, Male Luer Slip



WA PART- NUMBER	PRICE LIST DESCRIPTION
9209	Adapter, Single Male Luer Slip to Double Female Luer Slip, 1/EA
9226	Adapter, Single Female Luer Lock to Double Female Luer Slip , 1/EA
9229	Adapter, Double Female Threaded Connector to Double Female Luer Slip, 1/EA
9241	Adapter, Single Bayonet (Rectus) to Double Female Luer Slip , 1/EA
9281	Adapter, Double Female Threaded Connector to Double Female Luer Slip, 1/EA
9289	Adapter, Single 1/8" Hose Barb to Double Female Luer Slip, 1/EA
9304	Adapter, Single Male Luer Slip to Double Female Luer Slip , 1/EA
9357	Adapter, Double Female Subminiature to Double Female Luer Slip, 1/EA

Adapters for connecting cuffs with: 2 Tube, Female Luer Slip

WA PART- NUMBER	PRICE LIST DESCRIPTION
9074	Adapter, Single Female Threaded Connector to Double Male Slip Luer, 1/EA
9101	Adapter, Single 5/32" Hose Barb to Double Male Luer Slip, 1/EA
9205	Adapter, Double 1/8" Hose Barb to Double Male Luer Slip, 1/EA
9206	Adapter, Double Female Threaded Connector to Double Male Luer Slip, 1/EA
9207	Adapter, Double Female Subminiature to Double Male Luer Slip, 1/EA
9208	Adapter, Single Female Luer Lock to Double Male Slip Luer, 1/EA
9215	Adapter, Single Male Luer Lock to Double Male Slip Luer, 1/EA
9216	Adapter, Double Male Luer Slip to Double Male Luer Slip, 1/EA
9218	Adapter, Double Male Threaded Connector to Double Male Luer Slip, 1/EA
9220	Adapter, Single Male Luer Slip to Double Male Slip Luer, 1/EA
9224	Adapter, Single Female Subminiature to Double Male Slip Luer, 1/EA
9244	Adapter, Double Hose Barb & 10" Bare Tube to Double Male Luer Slip, 1/EA
9258	Adapter, Single Male Subminiature to Double Male Slip Luer, 1/EA
9325	Adapter, Double 1/8" Hose Barb & 10" Bare Tube to Double Male Luer Slip, 1/EA
9340	Adapter, Single 18" Hose Barb to Double Male Slip Luer, 1/EA
TL9203	Adapter, Double Mated Subminiature to Double Male Luer Slip, 1/EA
TL9204	Adapter, Single Bayonet (Rectus) to Double Male Slip Luer, 1/EA



Adapters for connecting cuffs with: 2 Tube, Female Subminiature



WA PART- NUMBER	PRICE LIST DESCRIPTION
9240	Adapter, Double Female Threaded Connector to Double Male Subminiature, 1/EA
9245	Adapter, Single Bayonet (Rectus) to Double Male Subminiature, 1/EA
9247	Adapter, Single Female Luer Lock to Double Male Subminiature, 1/EA
9268	Adapter, Single Female Threaded Connector to Double Male Subminiature, 1/EA
9305	Adapter, Single 5/32" Hose Barb to Double Male Subminiature, 1/EA

Adapters for connecting cuffs with: 2 Tube, Female Screw (Threaded) Connectors

WA PART- NUMBER	PRICE LIST DESCRIPTION
66044	Adapter, Single Female Luer Lock to Double Male Threaded Connector, 10/PK
66047	Adapter, Single Male Luer Lock to Double Male Threaded Connector, 10/PK
66048	Adapter, Single Female Threaded Connector to Double Male Threaded Connector, 10/PK
66050	Adapter, Single Bare Tube to Double Male Threaded Connector, 10/PK
66053	Adapter, Single Bayonet (Rectus) to Double Male Threaded Connector, 10/PK
9097	Adapter, Double Mated Subminiature to Double Male Threaded Connector, 1/EA
9230	Adapter, Single Bayonet (Rectus) to Double Male Threaded Connector, 1/EA
9231	Adapter, Single Female Luer Lock to Double Male Threaded Connector, 1/EA
9236	Adapter, Double Female Subminiature to Double Male Threaded Connector, 1/EA
9249	Adapter, Single 5/32" Hose Barb to Double Male Threaded Connector, 1/EA
9269	Adapter, Single Female Threaded Connector to Double Male Threaded Connector, 1/EA
9270	Adapter, Double 5/32" Hose Barb & Bare Tube to Double Male Threaded Connector, 1/EA
9291	Adapter, Single Female Subminiature to Double Male Threaded Connector, 1/EA
9292	Adapter, Double Female Luer & 5/32" Hose Barb to Double Male Threaded Connector, 1/EA
9293	Adapter, Single Male Luer Slip to Double Male Threaded Connector, 1/EA
9298	Adapter, Double Male Luer Slip to Double Male Threaded Connector, 1/EA
9302	Adapter, Single Male Subminiature to Double Male Threaded Connector, 1/EA
9334	Adapter, Double Female Luer to Double Male Threaded Connector, 1/EA



Adapters for connecting cuffs with: 2 Tube, Mated Subminiature









WA PART- NUMBER	PRICE LIST DESCRIPTION
9303	Adapter, Single Bayonet (Rectus) to Double Mated Subminiature Connectors, 1/EA
9308	Adapter, Single Female Threaded Connector to Double Mated Subminiature Connectors, 1/EA
9336	Adapter, Double Female Threaded Connector to Double Mated Subminiature Connectors, 1/EA
9339	Adapter, Double Female Subminiature to Double Mated Subminiature Connectors, 1/EA



Adapters for connecting cuffs to Manual Sphygmomanometers

WA PART- NUMBER	PRICE LIST DESCRIPTION
9228	Adapter, Double Female Luer & Bulb and Valve to Double Male Luer Slip, 1/EA
9242	Adapter, Single Male Luer Lock to Double Female Luer Slip & Bulb and Valve, 1/EA
9267	Adapter, Single Female Luer Slip to Double Female Bayonet (Rectus), Metal & Bulb and Valve, 1/EA
9272	Adapter, Single Extension Tubing to Double Male Threaded Connector & Bulb and Valve, 1/EA
9273	Adapter, Double Bulb and Valve & 10' Bare Tube to Double Male Luer Slip, 1/EA
9274	Adapter, Double Female Luer & Bulb and Valve to Double Male Threaded Connector, 1/EA
9294	Adapter, Double Bulb and Valve & 5' Coiled Tubing to Double Male Threaded Connector, 1/EA
9312	Adapter, Single 5' Coiled Tubing to Double Female Bayonet (Rectus), Metal, & Bulb and Valve, 1/EA
9313	Adapter, Single Female Luer Slip to Double Female Bayonet (Rectus) & Bulb and Valve, 1/EA
9317	Adapter, Double Male Luer Slip & Bulb and Valve to Double Mated Subminiature Connectors, 1/EA
9318	Adapter, Single Female Luer Slip to Double Male Slip Luer & Bulb and Valve, 1/EA
9320	Adapter, Single Female Luer Slip to Double Female Subminiature & Bulb and Valve, 1/EA
9330	Adapter, Single 5' Coiled Tubing to Double Female Bayonet (Rectus) & Bulb and Valve, 1/EA
9342	Adapter, Single 10" Bare Tube to Double Female Bayonet (Rectus), Metal & Bulb and Valve, 1/EA
TL9235	Adapter, Double Female Luer & Bulb and Valve to Double Male Subminiature, 1/EA
TL9271	Adapter, Single Female Luer Slip to Double Male Threaded Connector & Bulb and Valve, 1/EA

Adapters for connecting cuffs to Manual Sphygmomanometers





- 9242 No Illustration
- 9272 No Illustration
- 9273 No Illustration
- 9274 No Illustration
- 9317 No Illustration
- 9318 No Illustration
- 9320 No Illustration
- 9342 No Illustration



Adapters for Special Applications

WA PART-NUMBER	TRIMLINE DESCRIPTION
9114	ADAPTER,MONITOR,1-TUBE CABLE W/MALE LOCKING LUER CONNECTOR
9233	ADAPTER, DOUBLE TUBE W/ MALE SLIP LUERS, ONE FEMALE SLIP LUER AND ONE BARB
9296	ADAPTER, MANUAL SPHYG TO SGL TUBE BP CUFF
9301	ADAPTOR, MANUAL SPHYG TO SGL TUBE CUFF WITH BAYONET CNCTOR
9326	ADAPTER, MANUAL SPHYG W/ ONE TUBE, COILED TUBING, MALE SLIP LUER
9331	ADAPTER, 45" TUBING WITH MALE QUICK COUPLER CONNECTOR
9332	ADAPTER, 65" GREY TUBING WITH MALE QUICK COUPLER CONNECTOR
9354	ADAPTER,2-TUBE,2 MALE THREADED TO FEMALE LUER AND LENGTH OF TUB ING



9114 No Illustration
 9233 No Illustration
 9296 No Illustration
 9301 No Illustration
 9326 No Illustration
 9354 No Illustration