

Welch Allyn, Inc
4341 State Street Road
P.O. Box 220
Skaneateles Falls, NY 13153-0220
USA
Telephone: 315-685-4100



Re: Life-Cycle Testing on FlexiPort Reusable, Disposable, and EcoCuffs

To Whom It May Concern:

Welch Allyn does extensive lifecycle testing on our FlexiPort cuffs to ensure they will meet the demanding environment in which they are used. The cuff materials have been selected based on strength and patient comfort. To test the life of the cuff, Welch Allyn applies the cuffs to a mandrel to simulate a patient arm, inflates the cuff to 300 mmHg, deflates the pressure to 0 mmHg, and disengages the Velcro. This entire process is considered one cycle.

- The reusable FlexiPort cuffs are tested to over 30,000 cycles.
- The disposable FlexiPort cuffs are tested to over 500 cycles.
- The Ecocuffs are tested to over 100 cycles.

Please feel free to contact me if you have any additional questions.

Sincerely,

A handwritten signature in black ink that reads "Jamie Bartenstein".

Jamie Bartenstein
Global Category Specialist, Blood Pressure Cuffs
E: BartensteinJ@welchallyn.com
P: 315-685-2502