

Welch Allyn, Inc
4341 State Street Road
P.O. Box 220
Skaneateles Falls, NY 13153-0220
USA
Telephone: 315-685-4100



Re: FlexiPort Reusable Cuffs and Cleaning

To Whom It May Concern:

Welch Allyn has conducted extensive testing on the FlexiPort® reusable cuffs to determine the impact of using various cleaning solution and methods. Most cleaning solutions have minimal impact on the functionality of the cuffs. However, after periods of repetitive exposure, some solutions can cause cosmetic changes to the appearance of the device.

The cleaning solutions and methods identified in the instruction manual that ships with the product have been tested and validated. However, Welch Allyn realizes that some hospitals will use alternative solutions due to policies and availability. If the solution that your facility uses, specifically Caviwipes, is acceptable to your hospital and you are satisfied with the antimicrobial effectiveness of the product in your operating environment, then it is acceptable to use this product.

Sincerely,

A handwritten signature in black ink that reads "Jamie Bartenstein".

Jamie Bartenstein
US/C Acute Care Marketing, Blood Pressure Cuffs
E: Jamie.Bartenstein@welchallyn.com
P: 315-685-2502